

Sharing CGM data with Dexcom Clarity




Sharing your Dexcom CGM data with a clinic

Sharing data with a practice from the Dexcom CGM app

You can share data with your team using your Dexcom ONE+, G6 or G7 app.

Steps:


1. Open Dexcom ONE+ app > select 'connections' from the bottom menu
2. Select 'Clarity Clinic'
3. Enter clinic code: _____
4. Select 'confirm'
5. Confirm the clinic > Select 'confirm'
6. Select 'done' to return to the Home screen




Sharing CGM data with Dexcom Clarity




Sharing your Dexcom CGM data with a clinic




3. Enter clinic code:



4. Select 'continue'



5. Confirm the clinic
> Select 'confirm'



6. Select 'Done' to return Home

